



ADVAKIT Hernia

Ultimate strength and confidence for repair of abdominal wall hernia and deficiencies

ADVAKIT Hernia KIT
STERILISED SURGICAL KIT

AH01-01
MONOFILAMENT KNITTED POLYPROPYLENE UNDYED MESH
6 cm X 11 cm

ADVA1326
MONOFILAMENT POLYGLYCAPRONE 25 UNDYED
3-0 (2 metric)
70 cm

ADVA841
MONOFILAMENT POLYPROPYLENE BLUE
2-0 (3 metric)
70 cm

ADV0
MONOFILAMENT POLYPROPYLENE BLUE
0 (3.5 F)
70

Manufactured in India by
AMS
Advanced MedTech Solutions Pvt. Ltd. | Advanced MedTech Solutions Pvt. Ltd. | 21-22, 25-26 & 34-35, GPOC, Marjapur, tel: - 5600 08

AH01
STERILE
Net Quantity - 04 Nos.

25 mm 3/8 Circle CUTTING PRIME

30 mm 1/2 Circle ROUND BODIED

ADVA1326
3-0 (2 metric)
70 cm

ABSORBABLE SURGICAL SUTURE (SYNTHETIC)
ADVAMRYL
MONOFILAMENT POLYGLYCAPRONE 25 UNDYED
STERILISED SURGICAL NEEDLED SUTURE

ADVAMRYL
MONOFILAMENT POLYGLYCAPRONE 25 UNDYED
STERILISED SURGICAL NEEDLED SUTURE

ADVA841
2-0 (3 metric)
70 cm

NON ABSORBABLE SURGICAL SUTURE U.S.P.
ADVALENE
MONOFILAMENT POLYPROPYLENE BLUE
STERILISED SURGICAL NEEDLED SUTURE

ADVA842
0 (3.5 metric)
70 cm

NON ABSORBABLE SURGICAL SUTURE U.S.P.
ADVALENE
MONOFILAMENT POLYPROPYLENE BLUE
STERILISED SURGICAL NEEDLED SUTURE

ADVAMESH
MONOFILAMENT KNITTED POLYPROPYLENE
TRISIAL RECONSTRUCTION PROTHESIS
STERILISED SURGICAL KIT

ADVAPO611
6 cm X 11 cm

Kit Code	Brand Name	Code	Mesh Size (cm)	Suture		Needles	Mesh Per Box	Foil Per Box
				Size	Length (cm)			
AH01-01	ADVAMESH	ADVAP0611	6 X 11	-	-	-	01	03
	ADVAMRYL	ADVA1326	-	3-0	70	25 mm, 3/8 Circle, Cutting, PRIME		
	ADVALENE	ADVA841	-	2-0	70	30 mm, 1/2 Circle, Round Bodied		
	ADVALENE	ADVA842	-	0	70	30 mm, 1/2 Circle, Round Bodied		
AH01-02	ADVAMESH	ADVAP7515	7.5 X 11	-	-	-	01	03
	ADVAMRYL	ADVA1326	-	3-0	70	25 mm, 3/8 Circle, Cutting, PRIME		
	ADVALENE	ADVA841	-	2-0	70	30 mm, 1/2 Circle, Round Bodied		
	ADVALENE	ADVA842	-	0	70	30 mm, 1/2 Circle, Round Bodied		
AH02-01	ADVAMESH	ADVAP0611	6 X 11	-	-	-	01	03
	ADVACRYL	ADVA2317	-	2-0	90	30 mm, 1/2 Circle, Round Bodied		
	ADVALENE	ADVA841	-	2-0	70	30 mm, 1/2 Circle, Round Bodied		
	ADVALON	ADVA3336	-	2-0	70	45 mm, 3/8 Circle, Reverse Cutting		
AH02-02	ADVAMESH	ADVAP7515	7.5 X 11	-	-	-	01	03
	ADVACRYL	ADVA2317	-	2-0	90	30 mm, 1/2 Circle, Round Bodied		
	ADVALENE	ADVA841	-	2-0	70	30 mm, 1/2 Circle, Round Bodied		
	ADVALON	ADVA3336	-	2-0	70	45 mm, 3/8 Circle, Reverse Cutting		
AH02-03	ADVAMESH	ADVAP7615	7.6 X 11	-	-	-	01	03
	ADVALENE	ADVA844	-	2-0	90	30 mm, 1/2 Circle, Round Bodied		
	ADVACRYL	ADVA2317	-	2-0	90	30 mm, 1/2 Circle, Round Bodied		
	ADVALON	ADVA3336	-	2-0	70	45 mm, 3/8 Circle, Reverse Cutting		
AH02-04	ADVAMESH	ADVAP0611	6 X 11	-	-	-	01	03
	ADVALENE	ADVA843	-	1	70	40 mm, 1/2 Circle, Round Bodied, HEAVY		
	ADVALON	ADVA3328	-	3-0	70	26 mm, 3/8 Circle, Reverse Cutting		
AH02-05	ADVAMESH	ADVAP0611	6 X 11	-	-	-	01	03
	ADVACRYL	ADVA2317	-	2-0	90	30 mm, 1/2 Circle, Round Bodied		
	ADVALON	ADVA3336	-	2-0	70	45 mm, 3/8 Circle, Reverse Cutting		
AH02-06	ADVAMESH	ADVAP0611	6 X 11	-	-	-	01	03
	ADVACRYL	ADVA2347	-	1	90	40 mm, 1/2 Circle, Round Bodied		
	ADVALENE	ADVA840	-	1	70	30 mm, 1/2 Circle, Round Bodied, HEAVY		
	ADVALON	ADVA3336	-	2-0	70	45 mm, 3/8 Circle, Reverse Cutting		
AH03-01	ADVAMESH	ADVAP7515	7.5 X 15	-	-	-	01	03
	ADVACRYL	ADVA2317	-	2-0	90	30 mm, 1/2 Circle, Round Bodied		
	ADVAMRYL	ADVA1326	-	3-0	70	25 mm, 3/8 Circle, Cutting, PRIME		
	ADVALENE	ADVA841	-	2-0	70	30 mm, 1/2 Circle, Round Bodied		
AH03-02	ADVAMESH	ADVAP1010	10 X 10	-	-	-	01	03
	ADVACRYL	ADVA2317	-	2-0	90	30 mm, 1/2 Circle, Round Bodied		
	ADVAMRYL	ADVA1326	-	3-0	70	25 mm, 3/8 Circle, Cutting, PRIME		
	ADVALON	ADVA3336	-	2-0	70	45 mm, 3/8 Circle, Reverse Cutting		

