



## ADVAKIT Ortho

Supreme handling performance for complete ease and confidence in Orthopaedic surgeries




Kit Code	Brand Name	Code	Suture		Needles	Foil Per Kit	Kit Per Box
			Size	Length (cm)			
A001	ADVACRYL	ADVA2421	1	90	40 mm, 1/2 Circle, Reverse Cutting, HEAVY	01	06
	ADVALON	ADVA3336	2-0	70	45 mm, 3/8 Circle, Reverse Cutting	01	
A002	ADVACRYL	ADVA2421	1	90	40 mm, 1/2 Circle, Reverse Cutting, HEAVY	01	06
	ADVALON	ADVA3328	3-0	70	26 mm, 3/8 Circle, Reverse Cutting	01	
A003	ADVACRYL	ADVA2382	2-0	90	40 mm, 1/2 Circle, Reverse Cutting	01	06
	ADVALON	ADVA3336	2-0	70	45 mm, 3/8 Circle, Reverse Cutting	01	
A004	ADVACRYL	ADVA2382	2-0	90	40 mm, 1/2 Circle, Reverse Cutting	01	06
	ADVALON	ADVA3328	3-0	70	26 mm, 3/8 Circle, Reverse Cutting	01	
A005	ADVACRYL	ADVA2534	0	90	36 mm, 1/2 Circle, Reverse Cutting	01	06
	ADVALON	ADVA3337	0	70	45 mm, 3/8 Circle, Reverse Cutting	01	
A006	ADVACRYL	ADVA2382	2-0	90	40 mm, 1/2 Circle, Reverse Cutting	01	06
	ADVACRYL	ADVA2421	1	90	40 mm, 1/2 Circle, Reverse Cutting, HEAVY	01	

## ADVACRYL Polyglactin - 910 (Undyed & Violet) Braided Coated Synthetic Absorbable Surgical Suture



Knot pull tensile strength (KgF)



## ADVALON Polyamide Monofilament Synthetic Non-Absorbable Surgical Suture



Knot pull tensile strength (KgF)

