

ADVASTAP Reloads with Tri-Sure Technology.



**ADVASTAP ENDO RELOAD
(Purple)**



**ADVASTAP ENDO RELOAD
(Black)**



**ADVASTAP ENDO RELOAD
(Tan)**

Codes	ADVAER30PN, ADVAER45PN, ADVAER60PN	ADVAER45BN, ADVAER60BN	ADVAER30TN, ADVAER45TN, ADVAER60TN
Description	Manual stapler reload	Manual stapler reload	Manual stapler reload
Device Technology	Gripping surface technology	Gripping surface technology	Gripping surface technology
Sizes	30mm, 45mm, 60mm	45mm, 60mm	30mm, 45mm, 60mm
Tissue Thickness Description	Medium/Thick	Thick/Extra Thick	Thick/Vascular
Closed Staple Height	1.25-1.5-1.75	1.75-2.0-2.25	0.75-1.0-1.25
Open Staple Height	3.0-3.5-4.0mm	4.0-4.5-5.0mm	2.0-2.5-3.0mm
Staple Type	B-Form	B-Form	B-Form
No. of Reloads (Per Box)	6	6	6



**ADVASTAP ENDO RELOAD
(Purple)**



**ADVASTAP ENDO RELOAD
(Black)**



**ADVASTAP ENDO RELOAD
(Tan)**

Codes	ADVAER30P, ADVAER45P, ADVAER60P	ADVAER45B, ADVAER60B	ADVAER30T, ADVAER45T, ADVAER60T
Description	Manual stapler reload	Manual stapler reload	Manual stapler reload
Device Technology	Gripping surface technology	Gripping surface technology	Gripping surface technology
Sizes	30mm, 45mm, 60mm	45mm, 60mm	30mm, 45mm, 60mm
Tissue Thickness Description	Medium/Thick	Thick/Extra Thick	Thick/Vascular
Closed Staple Height	1.25-1.5-1.75	1.75-2.0-2.25	0.75-1.0-1.25
Open Staple Height	3.0-3.5-4.0mm	4.0-4.5-5.0mm	2.0-2.5-3.0mm
Staple Type	B-Form	B-Form	B-Form
No. of Reloads (Per Box)	6	6	6



ADVASTAP POWERED RELOAD (Purple)



ADVASTAP POWERED RELOAD (Black)



ADVASTAP POWERED RELOAD (Tan)

Codes	• ADVAPER30PN, ADVAPER45PN, ADVAPER60PN	• ADVAPER45BN, ADVAPER60BN	• ADVAPER30TN, ADVAPER45TN, ADVAPER60TN
Description	Powered stapler reload	Powered stapler reload	Powered stapler reload
Device Technology	Gripping surface technology	Gripping surface technology	Gripping surface technology
Sizes	30mm, 45mm, 60mm	45mm, 60mm	30mm, 45mm, 60mm
Tissue Thickness Description	Medium/Thick	Thick/Extra Thick	Thick/Vascular
Closed Staple Height	1.25-1.5-1.75	1.75-2.0-2.25	0.75-1.0-1.25
Open Staple Height	3.0-3.5-4.0mm	4.0-4.5-5.0mm	2.0-2.5-3.0mm
Staple Type	B-Form	B-Form	B-Form
No. of Reloads (Per Box)	• 6	• 6	• 6



ADVASTAP POWERED RELOAD (Purple)



ADVASTAP POWERED RELOAD (Black)



ADVASTAP POWERED RELOAD (Tan)

Codes	• ADVAPER30P, ADVAPER45P, ADVAPER60P	• ADVAPER45B, ADVAPER60B	• ADVAPER30T, ADVAPER45T, ADVAPER60T
Description	Powered stapler reload	Powered stapler reload	Powered stapler reload
Device Technology	Gripping surface technology	Gripping surface technology	Gripping surface technology
Sizes	30mm, 45mm, 60mm	45mm, 60mm	30mm, 45mm, 60mm
Tissue Thickness Description	Medium/Thick	Thick/Extra Thick	Thick/Vascular
Closed Staple Height	1.25-1.5-1.75	1.75-2.0-2.25	0.75-1.0-1.25
Open Staple Height	3.0-3.5-4.0mm	4.0-4.5-5.0mm	2.0-2.5-3.0mm
Staple Type	B-Form	B-Form	B-Form
No. of Reloads (Per Box)	• 6	• 6	• 6

